

格重品质·宜医无忧

Quality oriented, concern surgeons' concerns.

**Geyi** 格宜



## 双筋定位专利 钉成型更完美

THE PATENT DOUBLE STEMS DESIGN EXACTLY  
LOCATION AND BETTERTFORMING OF STAPLES

注册证号:浙械注准20152020239

## 一次性使用管型痔吻合器

Disposable Circular Stapler for Hemorrhoid

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

## 关键特征及利益

- 1.更多的吻合钉, 更好的止血效果。
- 2.磨砂调节螺母, 避免打滑, 手感更清晰。
- 3.垫刀圈采用科技革新新材料, 切断时提示声更清脆, 让操作更安心。
- 4.精确的最小吻合间隙, 避免过度压迫造成组织损伤。
- 5.器械经表面处理, 亮度适中, 不因发光而刺眼。



## Key Features & benefits

1. More staples, more effective for stopping bleeding.
2. Frosted nut with wing and the adjustment function, avoid slipping and feeling is more perfect.
3. Using new scientific material making the reminder sound more clear so that the operation is secure.
4. Minimal gap in anastomosis avoids tissue damage caused by over pressing.
5. Well surface treatment in body part, moderate intensity, not glare because of reflection.

## 适用范围

产品适用于III°-IV°痔的痔切闭术、痔上粘膜环切术或直肠粘膜内套叠手术。

## Application Scope:

Mainly Suitable for III°-IV°hemorrhoid and the rectal mucosa in the intussusception operation

## 规格型号 Model&Specification

型号规格 Model&Specification	抵钉座直径 (mm) Pin base diameter	环形刀直径 (mm) Annular cutter diameter	钉高 (mm) High nail	钉数 (mm) The number of nails
GYHS-32	32	23	4.2±0.2	32
GYHS-34	34	25		34

### 浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房  
住所: 浙江省杭州市滨江区长河街道江虹路333号  
研祥科技大厦A幢7层701室  
电话: 0086-571-86676309  
传真: 0086-571-86676306  
网址: www.zigeyi.com  
邮箱: trade@zigeyi.com 邮编: 310051

### ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town, Jiande City, 311606, Zhejiang Province, P.R. China  
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang Zone, Hangzhou, 310052, Zhejiang Province, P.R. China  
Tel: 0086-571-86676309 Fax: 0086-571-86676306  
Web: www.geyi-medical.com  
E-Mail: trade@zigeyi.com / geiyimedical@gmail.com  
PC: 310051

**Geyi 格宜**