

格重品质·宜医无忧
Quality oriented, concern surgeons' concerns.

Geyi 格宜



注册证号: 浙械注准20172091049

一次性使用包皮切割缝合器

Disposable Circumcision Staplers

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

产品特性

1. 整体成形的钉砧,确保钉成形完整。
2. 垫刀圈采用科技革新材料,切断时提示音更清脆,让操作更安心。
3. 环形刀与推钉片一体化结构,确保切割缝合的完整。
4. 更大直径的螺纹设计,确保切割和缝合时拉力。
5. 人体工程学防滑设计,手感更舒适。

适用范围

适用于临床包皮切割缝合手术。

Application Scope:

Suitable for foreskin suture surgery.

Product Features

1. The overall forming of the anvil to ensure that the staples forming integrity.
2. Special structure of the mat ring, effectively prompted the cutting suture to complete .
3. The combination of the ring knife and the push piece piece ensures the integrity of the cutting suture.
4. Larger diameter thread design ensures tension during cutting and stitching.
5. Anti-skid ergonomic design, makes it comfortable to use.



规格型号 Model&Specification

型号规格 Model&Specification	凹形座内径 Inner diameter of Concave	环形刀直径 diameter of Circular knife	钉数 (mm) The number of nails
GYCD-9	9	13.5	10
GYCD-11	11	16.5	12
GYCD-15	15	21.7	15
GYCD-21	20	27.6	18
GYCD-27	26	33.2	20
GYCD-33	31	37	24

浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房
住所: 浙江省杭州市滨江区长河街道江虹路333号
研祥科技大厦A幢7层701室
电话: 0086-571-86676309
传真: 0086-571-86676306
网址: www.zigeyi.com
邮箱: trade@zigeyi.com 邮编: 310051

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town, Jiande City, 311606, Zhejiang Province, P.R. China
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang Zone, Hangzhou, 310052, Zhejiang Province, P.R. China
Tel: 0086-571-86676309 Fax: 0086-571-86676306
Web: www.geyi-medical.com
E-Mail: trade@zigeyi.com / geyimedical@gmail.com
PC: 310051

Geyi 格宜