

格重品质·宜医无忧

Quality oriented, concern surgeons' concerns.

Geyi 格宜



一次性使用电凝手术器械

Disposable Laparoscopic Electrocoagulation surgical instruments

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

一次性使用双极电凝切割钳 Disposable Laparoscopic Bipolar Forceps



产品特性

1. 一次性使用, 无菌包装, 杜绝交叉感染;
2. 多种钳头选择, 适用不同的场合;
3. 手柄人体工学设计, 兼容大小手, 握持舒适轻盈灵巧, 操作方便。

Product Features

1. Single use, EO Sterilized package, effectively prevent cross infection;
2. Multiple forceps jaws for choices, meet the demand of various application;
3. Ergonomic handle design, applicable for different sizes of hands, easy to hold and operate.

| 型号 Model | 规格 Specification | 外径 Outer Diameter | 工作长度 Working Length | 张开角度 Open Angle |
|-----------|------------------|-------------------|---------------------|-----------------|
| GYEE-BC02 | φ5 X330mm | 5mm | 330mm | ≥30° |

一次性使用电凝切割剪刀 Disposable Laparoscopic Monopolar Scissor



产品特性

1. 一次性使用, 无菌包装, 杜绝交叉感染;
2. 剪刀头锋利, 剪切病患组织轻便;
3. 头部高位绝缘, 防止误损伤;
4. 手柄人体工学设计, 兼容大小手, 握持舒适轻盈灵巧, 操作方便。

Product Features

1. Single use, EO Sterilized package, effectively prevent cross infection;
2. Sharp-edged scissor, easy to cut human organs and tissues;
3. Insulating tip, free from mistaken injury;
4. Ergonomic handle design, applicable for different sizes of hands, easy to hold and operate.

| 型号 Model | 规格 Specification | 外径 Outer Diameter | 工作长度 Working Length | 张开角度 Open Angle |
|-----------|------------------|-------------------|---------------------|-----------------|
| GYEE-DS01 | φ5 X330mm | 5mm | 330mm | ≥40° |

一次性使用电凝切割分离钳 Disposable Laparoscopic Monopolar Dissector



产品特性

1. 一次性使用, 无菌包装, 杜绝交叉感染,
2. 多种钳头选择, 适用不同的场合;
3. 手柄人体工学设计, 兼容大小手, 握持舒适轻盈灵巧, 操作方便。

Product Features

1. Single use, EO Sterilized package, effectively prevent cross infection;
2. Multiple forceps jaws for choices, meet the demand of various application;
3. Ergonomic handle design, applicable for different sizes of hands, easy to hold and operate.

| 型号 Model | 规格 Specification | 外径 Outer Diameter | 工作长度 Working Length | 张开角度 Open Angle |
|-----------|------------------|-------------------|---------------------|-----------------|
| GYEE-DF01 | φ5 X330mm | 5mm | 330mm | α≥50° |

一次性使用电凝切割抓钳 Disposable Laparoscopic Monopolar Grasper



产品特性

1. 一次性使用, 无菌包装, 杜绝交叉感染,
2. 多种钳头选择, 适用不同的场合;
3. 手柄人体工学设计, 兼容大小手, 握持舒适轻盈灵巧, 操作方便。

Product Features

1. Single use, EO Sterilized package, effectively prevent cross infection;
2. Multiple forceps jaws for choices, meet the demand of various application;
3. Ergonomic handle design, applicable for different sizes of hands, easy to hold and operate.

| 型号 Model | 规格 Specification | 外径 Outer Diameter | 工作长度 Working Length | 张开角度 Open Angle |
|-----------|------------------|-------------------|---------------------|-----------------|
| GYEE-FT01 | φ5 X330mm | 5mm | 330mm | α≥50° |
| GYEE-FT02 | φ5 X330mm | 5mm | 330mm | α≥50° |

一次性使用电凝切割器 Disposable Laparoscopic Monopolar L-Hook



产品特性

1. 一次性使用, 无菌包装, 杜绝交叉感染,
2. 多种电勾头选择, 适用不同的场合;
3. 手柄执笔式设计, 兼容大小手, 握持舒适轻盈灵巧, 精确操作。

Product Features

1. Single use, EO Sterilized package, effectively prevent cross infection;
2. Multiple forceps jaws for choices, meet the demand of various application;
3. Ergonomic handle design, applicable for different sizes of hands, easy to hold and operate.

| 型号 Model | 规格 Specification | 外径 Outer Diameter | 工作长度 Working Length |
|-----------|------------------|-------------------|---------------------|
| GYEE-CF01 | φ5 X330mm | 5mm | 330mm |

浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房
住所: 浙江省杭州市滨江区长河街道江虹路333号
研祥科技大厦A幢7层701室
电话: 0086-571-86676309
传真: 0086-571-86676306
网址: www.zigeyi.com
邮箱: trade@zigeyi.com 邮编: 310051

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town, Jiande City, 311606, Zhejiang Province, P.R. China
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang Zone, Hangzhou, 310052, Zhejiang Province, P.R. China
Tel: 0086-571-86676309 Fax: 0086-571-86676306
Web: www.geyi-medical.com
E-Mail: trade@zigeyi.com / geiyimedical@gmail.com
PC: 310051

Geyi 格宜