

格重品质·宜医无忧

Quality oriented, concern surgeons' concerns.

Geyi 格宜



注册证号码: 浙械注准20192060039

一次性使用内窥镜标本取物袋

Disposable Specimen Retrieval Bag

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

产品特性

1. 袋体透明柔软,采用生物相容性良好的高强度TPU材料。
2. 张开装置采用医用NI-TI记忆合金材料,袋口张开充分,方便标本放入。
3. 取出时自动收紧袋口,避免溢出造成肿瘤种植或组织液渗透引发的炎症,降低组织损伤风险。
4. 操作简单,安全有效。
5. 有多种规格袋体可供选择,满足不同手术需求。

适用范围

供临床微创窥镜手术下收集人体组织标本/异物
取出体外用

Application Scope:

For the collection of human tissue specimens/foreign bodies
for clinical use in endoscopic surgery

Product Features

1. The bag is transparent and soft, and adopts high-strength TPU material with good biocompatibility.
2. The opening device is made of medical NI-TI memory alloy material, and the mouth of the bag is fully opened to facilitate the insertion of specimens.
3. Automatically tighten the bag when removing, to avoid the inflammation caused by tumor implantation or tissue fluid penetration, and reduce the risk of tissue damage.
4. Simple, safe and effective.
5. There are a variety of bag sizes to choose from to meet different surgical needs.



规格型号 Model&Specification

型号规格 Model&Specification	规格 Specifications	容量/ml Capacity/ml
GYSF-I型	35	35
	80	80
	130	130
	200	200
	300	300
GYSF-II型	50	50
	100	100
	150	150
	250	250
	300	300

型号规格 Model&Specification	规格 Specifications	容量/ml Capacity/ml
GYSF-III型	50	50
	100	100
	150	150
	250	250
	300	300
GYSF-IV型	50	50
	100	100
	150	150
	250	250
	300	300

浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房
住所: 浙江省杭州市滨江区长河街道江虹路333号
研祥科技大厦A幢7层701室
电话: 0086-571-86676309
传真: 0086-571-86676306
网址: www.zigeyi.com
邮箱: trade@zigeyi.com 邮编: 310051

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town,
Jiande City, 311606, Zhejiang Province, P.R. China
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang
Zone, Hangzhou, 310052, Zhejiang Province, P.R. China
Tel: 0086-571-86676309 Fax: 0086-571-86676306
Web: www.geyi-medical.com
E-Mail: trade@zigeyi.com / geyimedical@gmail.com
PC: 310051

Geyi 格宜