

格重品质·宜医无忧
Quality oriented, concern surgeons' concerns.

Geyi 格宜

注册证号:浙械注准20192020406

一次性使用套管穿刺器

Disposable Trocar



规格 Specification	配置内容 Assembly			
Φ3	1支Φ3穿刺套管+1支Φ3穿刺针	1pc 3mm cannula with stopcock, 1pc 3mm dilating tip obturator		
Φ5	1支Φ5穿刺套管+1支Φ5穿刺针	1pc 5mm cannula with stopcock, 1pc 5mm dilating tip obturator		
Φ10	1支Φ10穿刺套管+1支Φ10穿刺针	1pc 10mm cannula with stopcock, 1pc 10mm dilating tip obturator		
Φ12	1支Φ12穿刺套管+1支Φ12穿刺针	1pc 12mm cannula with stopcock, 1pc 12mm dilating tip obturator		
Φ15	1支Φ15穿刺套管+1支Φ15穿刺针	1pc 15mm cannula with stopcock, 1pc 15mm dilating tip obturator		
套装A/Kit A	Φ5、Φ10穿刺套管各2支+ Φ5、Φ10穿刺针各1支+气腹针+标本袋	2pcs 5mm cannula with stopcock, 2pcs 10mm cannula with stopcock,	1pc 5mm dilating tip obturator 1pc 10mm dilating tip obturator	1pc veress needle, 1pc specimen bag
套装B/Kit B	2支Φ5、1支Φ10穿刺套管+ Φ5、Φ10穿刺针各1支+气腹针+标本袋	2pcs 5mm cannula with stopcock, 1pc 10mm cannula with stopcock,	1pc 5mm dilating tip obturator 1pc 10mm dilating tip obturator	1pc veress needle, 1pc specimen bag
套装C/Kit C	3支Φ5、1支Φ10穿刺套管+ Φ5、Φ10穿刺针各1支+气腹针+标本袋	3pcs 5mm cannula with stopcock, 1pc 10mm cannula with stopcock,	1pc 5mm dilating tip obturator 1pc 10mm dilating tip obturator	1pc veress needle, 1pc specimen bag
套装E/Kit E	1支Φ5、2支Φ10穿刺套管+ Φ5、Φ10穿刺针各1支+气腹针+标本袋	1pc of 5mm cannula with stopcock, 2pcs of 10mm cannula with stopcock,	1pc of 5mm dilating tip obturator 1pc of 10mm dilating tip obturator	1pc of veress needle, 1pc specimen bag
套装F/Kit F	2支Φ5、2支Φ12穿刺套管+ Φ5、Φ12穿刺针各1支+气腹针+标本袋	2pcs 5mm cannula with stopcock, 2pcs 12mm cannula with stopcock,	1pc 5mm dilating tip obturator 1pc 12mm dilating tip obturator	1pc veress needle, 1pc specimen bag
套装K/Kit K	2支Φ5穿刺套管+1支Φ5穿刺针	2pcs 5mm cannula with stopcock, 1pc 5mm dilating tip obturator		
套装L/Kit L	2支Φ10穿刺套管+1支Φ10穿刺针	2pcs 10mm cannula with stopcock, 1pc 10mm dilating tip obturator		
套装AA/Kit AA	1支Φ5、2支Φ12穿刺套管+ Φ5、Φ12穿刺针各1支+气腹针+标本袋	1pc 5mm cannula with stopcock, 2pcs 12mm cannula with stopcock,	1pc 5mm dilating tip obturator 1pc 12mm dilating tip obturator	1pc veress needle, 1pc specimen bag

浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房
住所: 浙江省杭州市滨江区长河街道江虹路333号
研祥科技大厦A幢7层701室
电话: 0086-571-86676309
传真: 0086-571-86676306
网址: www.zigeyi.com
邮箱: trade@zigeyi.com 邮编: 310051

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town,
Jiande City, 311606, Zhejiang Province, P.R. China
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang
Zone, Hangzhou, 310052, Zhejiang Province, P.R. China
Tel: 0086-571-86676309 Fax: 0086-571-86676306
Web: www.geyi-medical.com
E-Mail: trade@zigeyi.com / geymedical@gmail.com
PC: 310051

Geyi 格宜

浙江格宜医疗器械有限公司
ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

产品特性

1. 无需转换器。
2. 方便5-12M各器械自由进出。
3. 套管螺纹设计, 腹壁损伤最小化。
4. 锋利薄刃, 减小创伤。
5. 一次性使用, 有效避免交叉感染。

Product Features

1. No need for an extra reducer.
2. Universal seal, suitable for the access of all general laparoscopic instruments.
3. Threaded cannula maximizes abdominal wall retention.
4. Sharp blade, reduce damage.
5. Single use, effectively prevent cross infection.



一次性使用套管穿刺器优点 Disposable Trocar Features

三级密封设计,
气密双保险。
Triple seals design,
airproof double insured

人体工程学设计,
操作舒适
Ergonomic design,
comfortable operation

螺纹套管设计,
腹壁损伤最小化。
Threaded cannula minimizes
abdominal wall retention

锋利薄刃穿刺头设计,
创伤小, 愈合快。
Sharp blade, ensures less
damage, quick recovery.

产品特性

1. 无需转换器。
2. 方便5-12M各器械自由进出。
3. 套管螺纹设计, 腹壁损伤最小化。
4. 钝性尖端针头设计, 创伤更小, 愈合更快。
5. 一次性使用, 有效避免交叉感染。

Product Features

1. No need for an extra reducer.
2. Universal seal, suitable for the access of all general laparoscopic instruments.
3. Threaded cannula minimizes abdominal wall retention.
4. Blunt tip design, ensures less damage, quick recovery.
5. Single use, effectively prevent cross infection.

