



Zhejiang Geyi Medical Instrument Co., Ltd

TEL: 0086-571-86676309

[www.geyi-medical.com](http://www.geyi-medical.com)

**Disposable Trocars  
(Single Port)**



**Single Port Surgery Access Devices**




# Introduction

Geyi single port trocar is a minimally invasive, multi-instrument access device designed to facilitate single incision laparoscopic surgery (SILS).

Certain laparoscopic procedures can be carried out via a single incision in the abdomen.

It contains a maximum of four instrument ports and two insufflation valves to support a range of four-instrument surgical procedures.



**Clinical  
application**



# Specification

Model	Specifications	Contents
GYTR-VII	Kit J	1pc 5mm cannula, 1pc 10mm cannula
	Kit K	2pcs 5mm cannula, 1pc 10mm cannula
	Kit L	2pcs 5mm cannula, 2pc 10mm cannula
	Kit Q	1pc 5mm cannula, 2pc 10mm cannula

# Components



# Features

1. Soft material, ensure a large floating angle for instruments.
2. In conjunction with wound protector, firmly connected, it's easy to assemble and disassemble.
3. Each cannula is designed with multi-layer sealing, ensures a perfect air proof.
4. The seal is made of high-strength medical-grade silicone material, which is resistant to laparoscopic instruments in and out.
5. Various choices in specifications, which meets the different demands in surgeries.
6. One surgical puncture hole, but multi cannula channels. Less traumatic, fast recovery. Best choice for scarless surgery.