

The disposable endoscopic stapler places triple staggered rows of titanium staples and simultaneously cuts and divides the tissue.

Universal handle is compatible with all reloading units.

Articulating reloads, can be rotated 45 degree left and right, provides more flexibility in narrow surgery space.

3 rows reloads provides equal height staples and stepped height staples.



Disposable Endoscopic Stapler		
Item Code	Shaft Length	
GYRC-60	60mm	short stapler
GYRC-100	160mm	standard stapler
GYRC-260	260mm	XL stapler





Articulating Reloads				
Item Code	Color	Suture Length (mm)	Anvil Height (mm)	Anvil Quantity
GYRCZ-B30L	white	31	2.5	42
GYRCZ-B30H	purple	31	3.5	42
GYRCZ-B45L	white	46	2.5	66
GYRCZ-B45H	purple	46	3.5	66
GYRCZ-B45M	black	46	4.8	66
GYRCZ-B60L	white	61	2.5	90
GYRCZ-B60H	purple	61	3.5	90
GYRCZ-B60M	black	61	4.8	90

* Reloads are compatible with Endo Gia universal stapler

White reloads
(equal height staples)

Purple & Black reloads
(stepped height staples)



Disposable Endoscopic Stapler