

格重品质·宜医无忧

Quality oriented, concern surgeons' concerns.

Geyi 格宜



一次性使用套管穿刺器(甲状腺专用型)

Disposable Trocar (Thyroid Surgery)

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

产品特性

1. 双层密封设计, 气密性可靠;
2. 加长套管, 满足甲状腺手术需求;
3. 套管深螺纹设计, 固定更可靠。



Product Features

1. Double seal design, reliable air tightness.
2. Longer cannula, meets the demands of thyroid surgery.
3. Threaded cannula, perfectly fix during surgery.



规格型号 Model&Specification

产品型号 Model& Specification	规格 Specifications	配置内容 Contents
GYTR-I型	Φ5JL	1支Φ5JL穿刺套管+1支Φ5JL穿刺针 1pc 5mmJL cannula with stopcock, 1pc 5mmJL dilating tip obturator
	Φ10JL	1支Φ10JL穿刺套管+1支Φ10JL穿刺针 1pc 10mmJL cannula with stopcock, 1pc 10mmJL dilating tip obturator
	套装H1 Kit H1	2支Φ5JL、2支Φ10JL穿刺套管+1支Φ5JL、 1支Φ10JL穿刺针+气腹针+标本袋 2pcs 5mmJL cannula with stopcock 2pcs 10mmJL cannula with stopcock, 1pc 5mmJL dilating tip obturator 1pc 10mmJL dilating tip obturator 1pc veress needle, 1pc specimen bag
	套装H2 Kit H2	2支Φ5JL、1支Φ10JL穿刺套管+1支Φ5JL、 1支Φ10JL穿刺针+气腹针+标本袋 2pcs 5mmJL cannula with stopcock 1pcs 10mmJL cannula with stopcock, 1pc 5mmJL dilating tip obturator 1pc 10mmJL dilating tip obturator 1pc veress needle, 1pc specimen bag
	套装H3 Kit H3	3支Φ5JL、1支Φ10JL穿刺套管+1支Φ5JL、 1支Φ10JL穿刺针+气腹针+标本袋 3pcs 5mmJL cannula with stopcock 1pcs 10mmJL cannula with stopcock, 1pc 5mmJL dilating tip obturator 1pc 10mmJL dilating tip obturator 1pc veress needle, 1pc specimen bag
	套装H4 Kit H4	1支Φ5JL、2支Φ10JL穿刺套管+1支Φ5JL、 1支Φ10JL穿刺针+气腹针+标本袋 1pcs 5mmJL cannula with stopcock 2pcs 10mmJL cannula with stopcock, 1pc 5mmJL dilating tip obturator 1pc 10mmJL dilating tip obturator 1pc veress needle, 1pc specimen bag

浙江格宜医疗器械有限公司

生产地址: 浙江省建德市下涯镇湖塘路5号1号、2号厂房
住所: 浙江省杭州市滨江区长河街道江虹路333号
研祥科技大厦A幢7层701室
电话: 0086-571-86676309
传真: 0086-571-86676306
网址: www.zigeyi.com
邮箱: trade@zigeyi.com 邮编: 310051

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory Address: No.1,2 Factory, No.5 hutang Road, Xiaya Town,
Jiande City, 311606, Zhejiang Province, P.R. China
Headquarters: A-701, No. 333, Jianghong Road, Changhe Street, Binjiang
Zone, Hangzhou, 310052, Zhejiang Province, P.R. China
Tel: 0086-571-86676309 Fax: 0086-571-86676306
Web: www.geyi-medical.com
E-Mail: trade@zigeyi.com / geyimedical@gmail.com
PC: 310051

Geyi 格宜